AKSO[®]

quality you can measure

INSTRUCTION MANUAL

AK135 Max

EC BENCH METER

CONTENTS

1 - SPECIFICATIONS	4
2 - ACCESSORIES	6
3 - OVERVIEW	7
FRONT VIEW	7
DISPLAY	8
4 - INSTRUCTIONS OF USE	9
TURN ON – TURN OFF	9
MEASUREMENT - Conductivity/Resistivity/TDS/Salinity	9
ADJUSTMENT - EC	10
5 - ADDITIONAL FUNCTIONS	12
DISPLAY LIGHTING	12
SET CONDUCTIVITY CELL CONSTANT	12
TDS FACTOR ADJUSTMENT	13
FACTORY RESET	13
6 - SETTINGS	14
SETUP MODE	14
P01 Con - Set conductivity cell constant (K)	14
P02 CAL - Visualization of the specific constant	14
P03 SLP - Set TDS factor	14
P04 Adj - Set the offset adjustment for temperature measurement	14
P05 Und - Select the temperature measurement unit	15
P06 bL - Enable/disable backlight	15
P07 Clr - Inactive	15
P08 rSt - Restore factory settings	15
7 - MAINTENANCE	16
CLEANING THE CONDUCTIVITY CELL	16
DISPOSAL OF BATTERIES AND ELETRONICS	16

1 - SPECIFICATIONS

EC	Measurement range:	0.0 to 19.99 $\mu\text{S/cm}$ 20.0 to 199.9 $\mu\text{S/cm}$ 200 to 1999 $\mu\text{S/cm}$ 2.00 to 19.99 mS/cm 20.0 to 199.9 mS/cm
	Resolution:	0.01 $\mu\text{S/cm}$ 0.1 $\mu\text{S/cm}$ 1 $\mu\text{S/cm}$ 0.01 mS/cm 0.1 mS/cm
	Accuracy (using the appropriate K constant for the range):	$\pm (1\%\text{FS} + 1 \text{ digit})$
TDS	Measurement Range:	0.00 to 19.99 ppm 20.0 to 199.9 ppm 200 to 1999 ppm 2.00 to 19.99 ppt 20.0 to 199.9 ppt
	Resolution:	0.01 ppm 0.1 ppm 1 ppm 0.01 ppt 0.1 ppt
	Accuracy:	$\pm (1\%\text{FS} + 1 \text{ digit})$
Salinity	Measurement Range:	0.0 to 100.0 ppt
	Resolution:	0.1 ppt
	Accuracy:	$\pm (1\%\text{FS} + 1 \text{ digit})$
Resistivity	Measurement Range:	0.0 to 100.0 $\text{M}\Omega$
	Resolution:	0.1 $\text{M}\Omega$
	Accuracy:	$\pm (1\%\text{FS} + 1 \text{ digit})$
Constant K	Configurable:	0.05 to 0.20 0.80 to 1.20 9.80 to 10.20
	Recommendation for conductivity cell and measurement range	
	K = 0.1	0.0 to 19.99 $\mu\text{S/cm}$
K = 0.1 and K = 1.0		20.0 to 199.9 $\mu\text{S/cm}$
K = 1.0		200 to 1999 $\mu\text{S/cm}$
K = 1.0		2.00 to 19.99 mS/cm
K = 10.0		20.0 to 199.9 mS/cm

1 - SPECIFICATIONS

TDS factor	Configurable:	0.30 to 1.00
Temperature compensation (Meter)	Automatic:	0 to 100°C
	Manual:	0 to 100°C
Temperature compensation (Cell)		0 to 80°C
Backlight:		On/Off
Reading freeze (HOLD):		Manual
Operating temperature of the conductivity cell:		0 to 50 °C
Operating humidity:		10 to 90 %UR (non-condensing)
Level of protection:		IP65
Power:		220VAC/12VDC
Dimensions (WxHxD):		216 x 190 x 59 mm
Weight:		950g
Additional functions:		
- Automatic calculation of the specific K constant; - Automatic factory calibration		

2 - ACCESSORIES

AK135 Max

Items included with the instrument:

- 1 conductivity cell with temperature sensor $K = 1.0$
- 1 storage and transport case
- 1 instruction manual
- 1 sample of $1413 \mu\text{S/cm}$
- 1 sample of 12.88 mS/cm

Items sold separately:

- Calibration solution $1413 \mu\text{S/cm}$ 250mL – AK4523
- Calibration solution 12.88 mS/cm 250mL – AK4524
- Conductivity cell with temperature sensor $K = 0.1$ – AK1021
- Conductivity cell with temperature sensor $K = 10.0$ – AK1022

K = 0.1

AK1021

K = 1.0

AK1020

K = 10.0


AK1022

Before use, examine the instrument and accompanying items carefully. If you detect any abnormalities, contact AKSO.

3 - OVERVIEW

FRONT VIEW


AK135 Max:

- 1 – Display
- 2 – On/Off/Set button
- 3 – Mode button
- 4 – Up button
- 5 – Down button
- 6 – CAL/ESC button
- 7 – HOLD/ENTER button

3 - OVERVIEW

DISPLAY LCD

AK135 Max:

- | | |
|---|---|
| 1 - Inactive | 19 - Manual temperature compensation |
| 2 - Active adjustment/calibration mode | 20 - Automatic temperature compensation |
| 3 - Conductivity cell constant | 21 - Settings menu description |
| 4 - Active reading freeze | 22 - Inactive |
| 5 - Inactive | 23 - Inactive |
| 6 - Battery replacement required | 24 - Measurement stabilization notification |
| 7 - MΩ measurement unit | 25 - Measurement value |
| 8 - Parts per million measurement unit | |
| 9 - Parts per thousand measurement unit | |
| 10 - μS/cm measurement unit | |
| 11 - mS/cm measurement unit | |
| 12 - Temperature measurement unit in °C | |
| 13 - Conductivity measurement mode | |
| 14 - Temperature measurement unit in °F | |
| 15 - Resistivity measurement mode | |
| 16 - Temperature measurement value | |
| 17 - TDS measurement mode | |
| 18 - Salinity measurement mode | |


4 – INSTRUCTIONS OF USE

POWER ON – POWER OFF

- To turn the instrument on or off, press the button  ;

MEASUREMENT – Conductivity / Resistivity / TDS / Salinity

- 1) Connect the conductivity cell to the instrument;
- 2) Remove the tip protection from the cell by rotating it counterclockwise if it is a cup (cells $K = 0.10$ and $K = 10.00$) or pulling it if it is the $K = 1.0$ cell;

- 3) Turn on the product with  button;
- 4) Rinse the cell in distilled water and remove excess water using a soft paper towel.
- 5) Immerse the cell in the sample to be analyzed and homogenize it.
- 6) Using the  button, select the desired measurement mode. For each press, the display will show:
 - **EC:** electrical conductivity
 - **RES:** resistivity
 - **TDS:** total dissolved solids
 - **SAL:** salinity
- 7) Wait for the reading to stabilize. The display will show 😊;

4 – INSTRUCTIONS OF USE

- 8) After the reading to stabilize, check the measurement value;
- 9) After completing the measurements, rinse the cell with distilled water and store preferably dry.

NOTES:

- If a temperature sensor for automatic temperature compensation (**ATC**) is not connected, the instrument will perform manual temperature compensation (**MTC**), which can be modified using the button  and  ;
- The instrument has the auto-calibration function whenever it is manually restarted after restoring the factory settings. Therefore, calibration with standard solutions should only be performed if a deviation beyond accuracy is detected.

ADJUSTMENT - EC

The calibration with the standard solution should be performed only with the solution closest to the operational measurement range in which conductivity is being measured. If the instrument is used for measurements close to 0 $\mu\text{S/cm}$, choose a standard solution close to this value.

The instrument does not have a fixed point, which means that the user has the freedom to calibrate with the solution they deem necessary in a practical and quick manner.

To make the adjustment, follow the instructions below:

- 1) Connect the conductivity cell to the instrument;
- 2) Remove the protective cap from the tip of the cell;
- 3) Turn on the instrument by pressing the  button;
- 4) Ensure that the measurement selected on the display is electrical conductivity by checking the **EC** icon on the display (if the instrument is set to another parameter, the conductivity measurement value will automatically be displayed when entering the adjustment mode);
- 5) Press the  button to enter the adjustment mode;
- 6) The indicator  will appear at the top of the screen signaling that the meter is in calibration mode. Note that the displayed value refers to the conductivity measurement without any adjustments;
- 7) Rinse the conductivity cell with distilled water and remove the excess water using a soft paper towel;

4 – INSTRUCTIONS OF USE

- 8) Prepare a portion of the solution to be used in a separate container. Ensure there is enough volume to cover the tip of the cell. **NEVER** insert the electrode into the original 250mL containers.
- 9) Immerse the conductivity cell in the standard solution, gently stirring to homogenize it;
- 10) Wait for the reading to stabilize. Then, using the buttons  and  adjust the measured value to the actual value of your solution. If necessary, press and hold the button to correct the value more quickly.
- 11) Once the reading has stabilized, press the  button to save the adjustment. The instrument will automatically return to measurement mode with the corrected value.

NOTE:

- *The adjustment of the conductivity parameter affects the other calculations performed by the meter.*

5 – ADDITIONAL FUNCTIONS

DISPLAY LIGHTING

To activate/deactivate the display lighting:


- 1) In measurement mode, press and hold the  button to enter the setup mode.
The display will show the indication P01;
- 2) Using the buttons  and  navigate to the parameter **P06 bL** and briefly press the  button to access it;
- 3) Using the buttons  and  select the option **On** to activate or **OFF** to deactivate. Press the  button to confirm.
- 4) Press the  button to return to measurement mode.

SET CONDUCTIVITY CELL CONSTANT

The conductivity cell constant is directly related to the measurement range in which the equipment will be used. To obtain stable and accurate results (ensuring the product's precision), it is necessary to use the cell according to the table below:

K = 0.1	0.0 a 19.99 $\mu\text{S/cm}$
K = 0.1 and K = 1.0	20.0 a 199.9 $\mu\text{S/cm}$
K = 1.0	200 a 1999 $\mu\text{S/cm}$
K = 1.0	2.00 a 19.99 mS/cm
K = 10.0	20.0 a 199.9 mS/cm

The exact value of the cell constant is usually displayed on the cable, as shown in the figure:

- 1) In measurement mode, press and hold the  button to enter the setup mode.
The display will show the indication P01;
- 2) Using the buttons  and  navigate to the parameter **P01 Con** and briefly press the  button to access it;
- 3) Using the buttons  and , increment or decrement the value according to your conductivity cell. Then, press the  button to confirm.
- 4) The product will return to configuration mode. Check if the option saved in memory is desired one.

5 – ADDITIONAL FUNCTIONS

TDS FACTOR ADJUSTMENT

The EC Max has the function of estimating the total dissolved solids in the solution by converting conductivity using a known constant. The TDS factor depends on the sample's reference salt and can be adjusted as needed by following the instructions below:

- 1) In measurement mode, press and hold the  button to enter the setup mode.
The display will show the indication P01;
- 2) Using the buttons  and , navigate to the parameter **P03 tdS** and briefly press the  button to access it;
- 3) Using the buttons  and , increment or decrement the value according to your reference salt. Then, press the  button to confirm.
- 4) The product will return to configuration mode. Check if the option saved in memory is the desired one.
- 5) Press the  button to return to measurement mode.

FACTORY RESET

Whenever an adjustment is made incorrectly or a complete restoration of the equipment's firmware is required, follow the instructions below:


- 1) In measurement mode, press and hold the  button to access the configuration mode. The display will show the indication P0;
- 2) Using the buttons  and , navigate to the parameter **P08 rSt** and briefly press the button  to access it;
- 3) "Yes" will be displayed on the screen. To confirm the reset, press the  button.
- 4) After the product automatically returns to the measurement screen, turn the instrument off and then on again.

6 - SETTINGS

SETUP MODE

- To access the instrument's configurable parameter menu, in measurement mode, press and hold the button . The display will show the indication P01;
- To navigate between parameters, use the buttons  e ;
- To access a parameter, press the button ;
- To toggle between the parameter configuration options use the buttons  and ;
- To confirm a setting, press the button . The instrument will return to the parameter selection;
- To go back to a previous level, press the button .

The figure below shows how the parameters are displayed:

Available parameters for configuration:

P01 Con – Set conductivity cell constant (K)

→ **0.10** (0.05 to 0.20)

→ **1.00** (0.80 to 1.20)

→ **10.00** (9.80 to 10.20)

P02 CAL – Visualization of the specific constant

This menu is dedicated to show the constant that the instrument is using in its measurements after user calibration.

P03 tdS – Set TDS factor

→ **0.30 to 1.00**

P04 Adj – Set the offset adjustment for temperature measurement

→ **± 5.0°C** (from the measured value)

NOTE:

- Ensure that the temperature sensor is connected to perform the temperature offset adjustment.

6 - SETTINGS

P05 Und – Select the temperature measurement unit

→ °C (Celsius)

→ °F (Fahrenheit)

P06 bL – Enable/disable backlight

→ On (enable)

→ OFF (disable)


P07 CLr – Inactive

→

P08 rSt – Restore factory settings

→ YES (perform o reset)

The figure below shows how the options for the accessed parameter are displayed:

1- Defined option

2- Another available option

3- Accessed parameter

Factory Settings of the Instrument

P01 Con	- Set Conductivity Cell Constant	1.00
P02 CAL	- Visualization of the Specific Constant	***
P03 tdS	- Set TDS factor	0.80
P04 AdJ	- Set Temperature Offset Adjustment	0.0
P05 Und	- Select Temperature Unit	°C
P06 bL	- Enable/disable backlight	OFF
P07 CLr	- Inactive	***
P08 rSt	- Reset	***

7 - MAINTENANCE

CLEANING THE CONDUCTIVITY CELL

To ensure the quality of measurements and the durability of the conductivity cell, make sure the tip of the cell is properly cleaned. For this, **NEVER** use:

- 1) Abrasive materials such as sandpaper, brushes;
- 2) Immersion or spraying of chemical solutions with extreme pH levels;
- 3) Aqueous solutions with temperatures higher than those specified in this document.

DISPOSAL OF BATTERIES AND ELECTRONICS

This product contains batteries and electronic components. Do not dispose of them with other common household waste. Deliver them to the appropriate collection center in accordance with local guidelines.

Important: the correct disposal of electronics and batteries prevents negative effects on the environment and for human health!

For more information about disposal services and locations, contact your local municipal office.

Akso guarantees its instruments against manufacturing defects with the following coverage: 2 years for bench meters, portable and pocket testers and 6 months for electrodes/sensors (unless otherwise indicated).

The warranty period starts from the original date of purchase by the customer. The warranty is only valid if the product is used under normal conditions and in accordance with its operating limits and recommended maintenance procedures.

WARRANTY

2

YEARS

This instrument comes with a 2-year* warranty against manufacturing defects and a 6-month* warranty for the sensor/electrode/probe.

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