
AKSO®

quality you can measure

INSTRUCTION MANUAL

AK90

POCKET PH METER

1. SPECIFICATIONS

	pH	Temperature
- Measuring range:	0.0 to 14.0	0.0 to 50.0°C
- Resolution:	0.1	0.1°C
- Accuracy:	± 0.1	± 1°C


- **Adjustment points:** pH 4.01, 7.00 and 10.01
- **Automatic temperature compensation (ATC):**
0 to 50 °C
- **IP rating:** IP54
- **Automatic shutdown:** 10 minutes
- **Operating temperature:** 0 to 50°C
- **Operating humidity:** 5 to 90%RH (non-condensing)
- **Power:** 6Vdc (4x LR44)
- **Dimensions (WxHxD):** 38 x 188 x 38mm
- **Weight:** 85g

Items sold separately:

- Electrode for AK90 AK1050
- Standard solution pH 4.01 AK4504
- Standard solution pH 7.00 AK4507
- Standard solution pH 10.01 AK4510
- KCl 3M storage solution AK4530
- Pepsin cleaning solution AK4531

IMPORTANTE: Before use, verify the AK90 carefully. If you detect any abnormalities, please contact AKSO.

2. OVERVIEW

3. INSTRUCTIONS OF USE

ON - OFF

To turn the AK90 on/off, briefly press the button **On/Off**.

ELECTRODE ACTIVATION

Before using you AK90 for the first time or after storing it dry (without KCL solution), activate the electrode according to the following procedure:

- 1) Remove the electrode protection cup;
- 2) Rinse the electrode with distilled water in abundance;
- 3) Fill approximately half of the electrode protection cup with 3 Molar potassium chloride solution (KCl-3M);
- 4) Insert the electrode into its protective cup, closing it tightly to prevent the KCl-3M solution from leaking;
- 5) Wait at least 30 minutes;
- 6) Make the adjustments to the AK90 measurement.

Your instrument will be ready for measurements.

MEASUREMENT ADJUSTMENT (CALIBRATION)

The AK90 can be adjusted to the following values:

2 points

- pH 7.00 and 4.01
- pH 7.00 and 10.01

3 points

- pH 7.00, 4.01 and 10.01

To make the adjustment, follow the procedure below:

- 1) Separate approximately 50mL of each pH buffer solution that will be used in the adjustment/calibration, pH 7.00, pH 4.01 and/or pH 10.01;
- 2) Remove the electrode protection cup;
- 3) Turn on the AK90 by pressing the **On/Off** button;
- 4) Rinse the electrode with distilled or deionized water in abundance and remove excess water using a soft paper towel;
- 5) Immerse the electrode in the pH 7.00 buffer solution, agitating it gently to homogenize it;
- 6) After the reading stabilizes, keep pressing the button **CAL** until the message **CAL** appears on the display. Release the button;
- 7) Wait for the message **SA** to appear on the display followed by **End**, signaling that the adjustment has been saved and completed.
- 8) To adjust the instrument at points pH 10.01 and/or pH 4.01, repeat the steps 4 to 7, using the corresponding buffer solutions.

NOTE:

- *After completing the adjustment procedure, discard the solutions used in running water .*

NEVER REUSE SOLUTIONS!

- *Adjust/calibrate the AK90 periodically or before measurements.*

MEASUREMENT

- 1) Remove the electrode protection cup;
- 2) Turn on the AK90 by pressing the **On/Off** button;

- 3) Rinse the electrode with distilled or deionized water in abundance and remove excess water using a soft paper towel;
- 4) Immerse the electrode in the sample solution, agitating it gently to homogenize it;
- 5) After the reading stabilizes, check the reading values for the sample pH and temperature on the display;
- 6) For a new measurement, repeat steps 3 to 5.

NOTE:

- **When** immerse the instrument into the sample, **NEVER** let the sample go above the protective cup seal. Submersion above this limit can irreversibly damage your AK90.
- At the end of measurements, always rinse electrode in distilled or deionized water in abundance and store it in the protective cup with clean KCl-3M.

4. ADDITIONAL FUNCTIONS

HOLD

To activate/deactivate Reading Freeze, press the **HOLD** button. **HOLD** will appear on the display while the function is active.

TEMPERATURE UNIT SELECTION - °C / °F

- 1) With the instrument turned off, keep pressing the **CAL** button and turn on the AK90. Release the **CAL** button;
- 2) Press the **CAL** button to switch between temperature measurements units °C or °F;
- 3) Press the **HOLD** button to confirm the selected unit showed on the display.

5. MAINTENANCE

RESTORATION OF FACTORY SETTINGS

- 1) With the AK90 connected in measurement mode, press and hold the **HOLD** button until the message **CLr** appears on the display. Release the button;
- 2) Perform the pH measurement adjustment.

ELECTRODE CLEANING

To ensure the quality of measurements and the durability of the electrode, perform periodic cleaning:

- 1) Prepare a container with electrode cleaning solution (pepsin);
- 2) Submerge the tip of the electrode in the solution for 20 to 25 minutes;
- 3) After this period, remove the electrode from the solution and rinse it thoroughly with distilled/deionized water;
- 4) Remove excess water using soft paper towels;
- 5) Store the electrode in a potassium chloride (KCl 3M) reservoir for 1 hour;
- 6) After resting, perform the instrument adjustment.

NOTES:

- *NEVER leave the electrode in the cleaning solution for more than 30 minutes.*
- *NEVER use abrasives (brushes, rough paper, etc.) to clean the electrode.*
- *AKSO offers a ready-to-use electrode cleaning solution. Contact AKSO for more information.*

ELECTRODE REPLACEMENT

- 1) Remove the electrode protection cup;
- 2) Rotate the attachment ring (on the top of the electrode) counterclockwise and remove it;
- 3) Detach the used electrode from the AK90;
- 4) Carefully install the new electrode for AK90, paying attention to its correct placement and connection;
- 5) Reattach the electrode attachment ring tightening it firmly (without applying force).

BATTERY REPLACEMENT

When the display shows "---" or the characters begin to blink, replace the used batteries with new ones following the procedure below:

- 1) Pull the battery compartment cover (located on the top of the AK90) and remove it;
- 2) Remove the used batteries from the compartment;
- 3) Install 4 new batteries, ensuring proper polarity;
- 4) Replace the compartment cover by pressing it firmly (without forcing it).

DISPOSAL OF BATTERIES AND ELECTRONICS

This product contains battery and electronic components. Do not dispose of them with other common household waste. Deliver them to the appropriate collection center in accordance with local guidelines. **Important:** The correct disposal of electronics and batteries prevents negative effects on the environment and on human health!

WARRANTY

2
YEARS


This instrument comes with a 2-year* warranty against manufacturing defects and a 6-month* warranty for the sensor/electrode/probe.

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