# AKSO<sup>®</sup>

quality you can measure

INSTRUCTION MANUAL

## AK131 Max

BENCH PH METER WITH GLP

# CONTENTS

<b>1 - SPECIFICATIONS</b> .....	<b>4</b>
<b>2 - ACCESSORIES</b> .....	<b>5</b>
<b>3 - OVERVIEW</b> .....	<b>6</b>
FRONT VIEW .....	6
DISPLAY LCD .....	7
<b>4 - INSTRUCTIONS OF USE</b> .....	<b>8</b>
POWER ON - POWER OFF .....	8
Measurement pH .....	8
Measurement ORP .....	9
<b>5 - CALIBRATION</b> .....	<b>10</b>
Adjustment pH .....	10
Adjustment ORP .....	11
<b>6 - ADDITIONAL FUNCTIONS</b> .....	<b>12</b>
MEASUREMENT RECORDS .....	12
Save records to memory .....	12
View records in memory .....	13
Delete all records in memory .....	13
COMMUNICATION WITH THE PRINTER .....	14
COMMUNICATION WITH THE PC .....	14
DISPLAY LIGHTING .....	15
CALIBRATION AT CUSTOM POINTS .....	15
SETTING THE DATE AND TIME .....	16
GLP FUNCTION .....	17
<b>7 - SETTINGS</b> .....	<b>18</b>
SETUP MODE .....	18
P01 buF - Select pH calibration standard .....	18
P02 CUSt - Customize pH set points .....	18
P03 SLP - View the slope values from the last calibration .....	18
P04 AdJ - Set the offset adjustment for temperature measurement .....	18
P05 Und - Select the temperature measurement unit .....	18
P06 Aut - Automatic shutdown options .....	18
P07 Add - Enable/disable automatic records .....	18
P08 dAtE - View calibration expiration date .....	19
P09 Auto - Enable/disable automatic freezing of reading .....	19
P10 GLP - View information on good laboratory practices .....	19
P11 rTc - Date and time setting .....	19
P12 bL - Enable/disable backlight .....	19
P13 Clr - Delete all records from memory .....	19
P14 rSt - Restore factory settings .....	19
<b>8 - MAINTENANCE</b> .....	<b>20</b>
ELECTRODE CLEANING - pH/ORP .....	20
DISPOSAL OF BATTERIES AND ELECTRONICS .....	21

## 1 - SPECIFICATIONS

pH	<b>Measurement range:</b>	-2.000 to 19.999 pH / - 2.00 to 19.99 pH
	<b>Resolution:</b>	0.001 pH / 0.01 pH
	<b>Accuracy:</b>	$\pm 0.01$ pH / $\pm 0.002$ pH
ORP	<b>Measurement range:</b>	-1999.9 to 1999.9 mV
	<b>Resolution:</b>	0.1mV
	<b>Accuracy:</b>	$\pm 1$ mV
Temperature	<b>Measurement range:</b>	-5 to 100°C
	<b>Resolution:</b>	0.1°C
	<b>Accuracy:</b>	$\pm 0.5$ °C
Calibration (pH)	<b>USA Standards:</b>	1.68, 4.01, 7.00, 10.01, 12.45 pH
	<b>NIST Standards:</b>	1.68, 4.00, 6.86, 9.18, 12.45 pH
Temperature compensation	<b>Automatic:</b>	-5 to 100°C
	<b>Manual:</b>	-5 to 100°C
Memory:		1000 records
Communication with the PC:		Type A USB
Communication with the printer:		RS232
Display Lighting (Backlight):		On/Off
Reading freeze (HOLD):		Manual/Automatic
Operating temperature:		0 to 50 °C
Operating humidity:		10 to 90 %UR (non-condensing)
Level of protection:		IP65
Power:		220VAC/12VDC
Low battery indicator:		Visual
Dimensions (WxHxD):		216 x 190 x 59 mm
Weight:		950g

## 2 - ACCESSORIES

---

### **AK131 Max:**

#### **Items included with the instrument:**

- 1 double junction glass pH electrode – AK1005
- 1 Storage and transport case
- 1 Instruction manual
- 1 sample of pH 4 solution
- 1 sample of pH 7 solution
- 1 sample of pH 10 solution
- 1 sample of KCl solution
- 1 power supply 220VAC/12VDC – AK1075
- 1 Electrode holder for bench meters - AK1070
- 1 Temperature probe - AK1040
- 1 USB Cable - AK1090

#### **Items sold separately:**


- pH 4.01 Standard Solution (250 mL bottle) – AK4504
- pH 7.00 Standard Solution (250 mL bottle) – AK4507
- pH 10.01 Standard Solution (250 mL bottle) – AK4510
- 3M KCl Storage Solution (250 mL bottle) – AK4530
- ORP electrode BNC connector – AK1010

Before use, examine the instrument and accompanying items carefully. If you detect any abnormalities, contact AKSO.

### 3 - OVERVIEW

---

#### FRONT VIEW


#### AK131 Max:

- 1 – Display
- 2 – Cal/Esc button
- 3 – GLP/Enter button
- 4 – Down button
- 5 – Up button
- 6 – On/Off/Set button
- 7 – Mode/Print button

### 3 - OVERVIEW

#### DISPLAY LCD

#### AK131 Max:

- |                                         |                                           |
|-----------------------------------------|-------------------------------------------|
| 1 - GLP indication active               | 12 - Temperature measurement unit in °F   |
| 2 - Adjustment/calibration mode active  | 13 - Manual temperature compensation      |
| 3 - USB connectivity indication         | 14 - Automatic temperature compensation   |
| 4 - Reading freeze active               | 15 - Active log indication                |
| 5 - Slope display indication            | 16 - Active custom point indication       |
| 6 - Measurement value                   | 17 - Expired calibration indication       |
| 7 - ORP measurement unit                | 18 - Full memory indication               |
| 8 - pH measurement unit                 | 19 - Print indication                     |
| 9 - Slope calculation unit              | 20 - Settings access indication           |
| 10 - Temperature measurement value      | 21 - Calibrated set points                |
| 11 - Temperature measurement unit in °C | 22 - Measurement stabilization indication |

## 4 – INSTRUCTIONS OF USE

---

### POWER ON – POWER OFF

Always connect the power supply to the instrument first and then to the **220V** electrical outlet. **The power supply is not dual voltage.**

- 1) To turn the instrument on or off, press the button  ;
- 2) The following will be displayed on the screen, respectively:  
**Display self-test → Model/Version/Greeting → Day → Month → Year**
- 3) The product will start up on the measurement screen.

### Measurement - pH

- 1) Connect the pH electrode and temperature sensor/probe to the instrument;
- 2) Remove the KCl bottle from the electrode tip, turning it counterclockwise;
- 3) Turn on the instrument by pressing the button  ;
- 4) Select pH measurement by pressing the button  . The display will show **pH**;
- 5) Rinse the electrode (and temperature probe) with distilled/deionized water and remove excess water using a soft paper towel;
- 6) Immerse the electrode (and temperature probe) in the sample to be analyzed;
- 7) Wait for the reading to stabilize. The display will show 😊;
- 8) After the reading stabilizes, check the pH and temperature values;
- 9) After completing the measurements, rinse the electrode (and temperature probe) with distilled water and store the electrode in the storage bottle with potassium chloride solution (KCl – 3M).

#### **NOTE:**

- *If a temperature sensor/probe for automatic temperature compensation (ATC) is not connected, the instrument will perform manual temperature compensation (MTC), which can be modified using the button  and .*

## 4 – INSTRUCTIONS OF USE

---

### Measurement - ORP

- 1) Connect the ORP electrode to the instrument;
- 2) Remove the KCl bottle from the electrode tip;
- 3) Turn on the instrument by pressing the button  ;
- 4) Select the ORP measurement mode by pressing the button  . The display will show the **mV** indicator;
- 5) Rinse the electrode with distilled/deionized water and remove excess water using a soft paper towel;
- 6) Immerse the electrode in the sample to be analyzed;
- 7) Wait for the reading to stabilize. The display will show 😊;
- 8) After the reading stabilizes, check the ORP value;
- 9) After completing the measurements, rinse the electrode with distilled water and store the electrode in the storage bottle with potassium chloride solution (KCl – 3M).

## 5 - CALIBRATION


---

### Adjustment - pH

Adjust the instrument's pH measurement upon receiving it and at least once a week. The instrument can be adjusted to 1, 2, 3, 4, or 5 fixed or customizable pH points.

**USA fixed points:** 1.68, 4.01, 7.00, 10.01, 12.45 pH (all customizable)

**NST fixed points:** 1.68, 4.01, 7.00, 10.01, 12.45 pH (all customizable)

- 1) Connect the pH electrode and temperature sensor/probe to the instrument;
- 2) Remove the KCl bottle from the electrode tip, turning it counterclockwise;
- 3) Turn on the instrument by pressing the button  ;
- 4) Ensure that the measurement selected on the display is pH by viewing the pH unit icon;
- 5) Press the button  to access the pH measurement adjustment mode;
- 6) The display will show an indicator  at the top of the screen, signaling that the meter is in calibration mode;
- 7) Rinse the electrode (and temperature probe) in distilled/deionized water and remove excess water using a soft paper towel;
- 8) Using a separate bottle, take a sufficient portion of the pH solutions, enough to cover the electrode tip. **NEVER** insert the electrode into the original vials of 250mL.
- 9) Immerse the electrode (and temperature probe) in the pH buffer solution, agitating it gently to homogenize it;
- 10) Wait for the reading to stabilize. The notification  will appear on the screen;
- 11) After the reading stabilizes, press the button  to save the adjustment; The instrument will flash the value of the standard solution followed by the slope and the message END;
  - To calibrate another pH point, repeat steps 7 to 10 using the corresponding buffer solution;
- 12) After completing the adjustment, the instrument will display the adjusted points at the bottom of the display;
- 13) To finish saving the adjusted points, press the button  ;
- 14) The message "Date Cal Day" will appear on the display, indicating that the day on which the calibration will expire must be selected;

## 5 - CALIBRATION

---

- 15) Use  and  buttons to select the appropriate day and press the  button to confirm;
- 16) The message "Date Cal Month" will appear on the display, indicating that the month in which the calibration will expire must be selected;
- 17) Use the  and  buttons to select the appropriate month and press the  button to confirm;
- 18) The message "Date Cal Year" will appear on the display, indicating that you must select the year in which the calibration will expire;
- 19) Use the  and  buttons to select the appropriate year and press the  button to confirm;
- 20) The instrument will return to the measurement screen;
- 21) After completing the work, wash the electrode (and temperature probe) with distilled/deionized water and store the electrode in the storage bottle with potassium chloride solution (KCl - 3M).

### NOTE:

- *If an adjustment expiration date prior to the date of use is selected, the measurement screen will automatically display an expired adjustment indication.*
- *The slope is calculated between points, so it will be displayed when there are at least two points fitted.*
- *The slope indicates the electrode's performance and should be between 80% and 115%. If the slope is below 80% or above 115% clean the electrode, refer to **7 - MAINTENANCE>ELECTRODE CLEANING – pH/ORP**.*
- *If the message **Err** appears on the display when trying to save an adjustment, check the pH solutions and the electrode and repeat the procedure;*

### Adjustment - ORP

The ORP electrode does not require adjustment. However, you can check the electrode's performance and the accuracy of the readings using an ORP standard solution.

## 6 - ADDITIONAL FUNCTIONS

### MEASUREMENT RECORDS

#### Save records to memory

The instrument has two recording modes: manual and automatic. To record **manually**, follow the instructions below:

- 1) After the reading has stabilized, press and hold the button  until **REC** appears on the display followed by the number of the record made.
- 2) The number of the record allocated in the memory will be displayed in place of the time.

#### NOTE:

- If the pH parameter is selected on the display, the pH and mV values will be stored in memory when recorded.
- If the ORP parameter is selected on the display (mV), only the mV value will be stored in memory.

To register **automatically**, follow the instructions below:

- 1) First, configure the interval and number of records via the settings menu. To do this, press and hold the  button;
- 2) In the menu, use the  and  buttons to navigate to the **P07 Add** option and briefly press the  button;
- 3) The message "YES" will appear on the display, press the  button;
- 4) The indication "Minute" will be displayed, use the buttons  and  to select the recording interval (1 to 120 minutes) and confirm with the  button;
- 5) The indication "Qty" will be displayed. Using the  and  buttons, select the number of records (1 to 200) and confirm with the  button.
- 6) After setting the recording interval and number of records, press and hold the  button on the measurement screen. The instrument will display the message REC and start automatic recording.
- 7) If you wish to cancel the recordings, briefly press the  button;
- 8) At the end of the recording period, the instrument will automatically return to measurement mode.

## 6 - ADDITIONAL FUNCTIONS

---

### NOTE:

- During the automatic recording period, the instrument will lock the use of all buttons except  and  ;
- If the pH parameter is selected on the display, the pH and mV values will be stored in memory when recorded.
- If the ORP parameter is selected on the display (mV), only the mV value will be stored in memory.

### View records in memory

To view measurement records saved in the instrument's memory:

- 1) In measurement mode, press and hold the  button until **REC** appears on the display followed by the number of the last recording made.
- 2) To navigate between records in memory, use the  and  buttons.
- 3) To switch between the pH and mV values of the same record, briefly press the  button.
- 4) To return to measurement mode, press the  button.

### Delete all records in memory

To delete all records stored in the instrument's memory:

- 1) In measurement mode, press and hold  button to access the configuration mode. P0 will appear on the display.
- 2) Using the  and  buttons, navigate to parameter **P13 CLr** and briefly press the  button to access it;
- 3) The message YES will be displayed. Press the  button to confirm. All records will be deleted and the instrument will return to the measurement screen.

## 6 - ADDITIONAL FUNCTIONS

### COMMUNICATION WITH THE PRINTER

To print a product registration, follow the instructions below:

1) You must use a compatible printer with the following characteristics:

Communication: Rs232

BaudRate: 9600 bps

Dimensions: 58 mm thermal

Cable: Serial DB9 M x DB9 F

2) Connect the serial cable to the printer and the instrument (printer connector);

3) Access the memory log display

(see **6 – ADDITIONAL FUNCTIONS>View memory logs**) and, on the desired log,

press and hold the  button to print, as indicated by the  icon on the screen.

### COMMUNICATION WITH THE PC

1) With the device turned **OFF**, connect the USB cable supplied with the product to the computer and the device (USB connector).

2) The device will turn on automatically and display the message “PC” and the USB indicator at the top of the display;

3) On the computer, open the device folder and, if there are records in the device's memory, you will see the file “memory.csv”.

4) The file will contain the records organized as follows:

Date	Time	pH	mV	Temp	Und	Calibrated Points						
08/09/2025	10:08:50	7.00	0.00	25.00	°C	0	0 4.01	0 7.00	0 10.01	0		
08/09/2025	10:09:50	4.01	170.86	25.00	°C	0	0 4.01	0 7.00	0 10.01	0		
08/09/2025	10:10:51	10.01	-172.0	25.00	°C	0	0 4.01	0 7.00	0 10.01	0		
08/09/2025	10:14:01	5.35	94.29	25.00	°C	0	0 4.01	0 7.00	0 10.01	0		
08/09/2025	10:14:24	9.18	-124.57	25.00	°C	0	0 4.01	0 7.00	0 10.01	0		

#### NOTE:

- Communication will only be successful if the device is turned off before connecting the cable, otherwise the message “Only OFF” will appear on the display indicating that the product must be turned off;
- During connection, it is not possible to take measurements or change any product settings.

## 6 - ADDITIONAL FUNCTIONS

---

### DISPLAY LIGHTING

To activate/deactivate the display lighting:

- 1) In measurement mode, press and hold the  button to enter the setup mode.  
The display will show the indication P01;
- 2) Using the buttons  and  navigate to the parameter **P12 bL** and briefly press the  button to access it;
- 3) Using the buttons  and  select the option **On** to activate or **OFF** to deactivate. Press the  button to confirm.
- 4) Press the  button to return to measurement mode.

### CALIBRATION AT CUSTOM POINTS

To calibrate at custom points, use the **P02 CUST** menu following the guidelines below:

- 1) In measurement mode, press and hold the  button to access the configuration mode. P0 will appear on the display;
- 2) Using the  and  buttons, navigate to parameter **P02 CUST** and briefly press the  button to access it;
- 3) Using the  and  buttons, select the set point you want to customize and briefly press the  button;
- 4) The value will flash on the display, indicating that it can be changed. Use the  and  buttons to enter the desired value and confirm with the  button.
- 5) Press the  button to return to measurement mode, and the next time you adjust the setting, the changed point will be used as the new default value.

## 6 - ADDITIONAL FUNCTIONS

---

### SETTING THE DATE AND TIME

To adjust the date and time, follow the instructions below:

- 1) In measurement mode, press and hold the  button to access the configuration mode. P0 will appear on the display;
- 2) Using the  and  buttons, navigate to parameter **P11 rtC** and briefly press the  button to access it;
- 3) Using the  and  buttons, select "date" if you want to set the date or "hour" to set the time, and press the  button;
- 4) If you select "hour," the time will flash on the display. Adjust using the  and  buttons and confirm with the  button;
- 5) Afterwards, the display will flash the indication "minute." Adjust using the  and  buttons and confirm with the  button.
- 6) If you select "date," the day will flash on the display. Adjust using the  and  buttons and confirm with the  button.
- 7) Afterwards, the "minute" will flash on the display. Adjust using the  and  buttons and confirm with the  button.
- 8) Next, repeat the process for the month and year. When finished, the product will return to the P11 menu. Press briefly  to return to the measurement screen.

## 6 - ADDITIONAL FUNCTIONS

---


### GLP FUNCTION

The GLP (Good Laboratory Practices) function stores information related to equipment calibration, including:

- Calibration expiration date;
- Electrode slope in each adjustment range;

This information can be checked in the **P10 GLP** menu and also by following the instructions below:

1) With the instrument turned on and on the measurement screen, press and hold the

2) The instrument will access and display the expiration date of the last valid adjustment.

To switch between date and slope information, use the  and  buttons;

**NOTE:**

- *If the product does not have a valid setting, the message "No Cal" will be displayed.*

## 7 - SETTINGS

---

### SETUP MODE

- To access the instrument's configurable parameter menu, in measurement mode, press and hold the button . The display will show the indication P0;
- To navigate between parameters, use the buttons  e ;
- To access a parameter, press the button ;
- To toggle between the parameter configuration options use the buttons  and ;
- To confirm a setting, press the button . The instrument will return to the parameter selection;
- To go back to a previous level, press the button .

Available parameters for configuration:

**P01 buF** – Select pH calibration standard

→ **USA** (1.68, 4.01, 7.00, 10.01 and 12.45 pH)

→ **NST** (1.68, 4.00, 6.86, 9.18 and 12.45 pH)

**P02 CUS<sub>t</sub>** – Customize pH set points

**P03 SLP** – View the slope values from the last calibration

**P04 Adj<sub>-</sub>** Set the offset adjustment for temperature measurement

→ **± 5.0°C** (from the measured value)

**NOTE:**

- *Ensure that the temperature sensor is connected to perform the temperature offset adjustment.*

**P05 Und** – Select the temperature measurement unit

→ **°C** (Celsius)

→ **°F** (Fahrenheit)

**P06 Aut** – Automatic shutdown options

→ **OFF** (Inactive automatic shutdown)

→ **5** (5 minutes)

→ **10** (10 minutes)

→ **30** (30 minutes)

→ **60** (60 minutes)

**P07 Add** – Enable/disable automatic records

→ **No** (enable)

→ **Yes** (disable)

## 7 - SETTINGS

**P08 dAtE** – View calibration expiration date

**P09 Auto**– Enable/disable automatic freezing of reading

→ **On** (enable)

→ **OFF** (disable)

**P10 GLP** – View information on good laboratory practices

**P11 rtC** – Date and time setting

**P12 bL** – Enable/disable backlight

→ **On** (enable)

→ **OFF** (disable)

**P13 Clr** – Delete all records from memory

→ **YES**

**P14 rSt**– Restore factory settings

→ **YES** (perform the reset)

### Factory Settings of the Instrument

<b>P01 buF</b>	- pH Calibration Standard	USA
<b>P02 CUSt</b>	- Customize pH set points	***
<b>P03 SLP</b>	- View Slope Values	***
<b>P04 AdJ</b>	- Set Temperature Offset Adjustment	0.0
<b>P05 Und</b>	- Select Temperature Unit	°C
<b>P06 Aut</b>	- Enable/disable automatic reading freeze	ON
<b>P07 Add</b>	- Enable/disable automatic records	NO
<b>P08 dAtE</b>	- View calibration expiration date	***
<b>P09 Auto</b>	- Enable/disable auto freezing of reading	ON
<b>P10 GLP</b>	- View information on GLP	***
<b>P11 rtC</b>	- Date and time setting	***
<b>P12 bL</b>	- Enable/Disable backlight	ON
<b>P13 Clr</b>	- Delete all records from memory	***
<b>P14 rSt</b>	- Reset	***

## 8 - MAINTENANCE

---

### ELECTRODE CLEANING - pH / ORP

To ensure measurement accuracy and prolong the life of the pH/ORP electrode, perform regular cleaning (at least every two weeks) or when the electrode slope falls below 90%:

- 1) Prepare a container with an electrode cleaning solution (5% pepsin in HCl 0.1M solution);
- 2) Immerse the electrode tip in the solution for 20 to 25 minutes;
- 3) After this period, remove the electrode from the solution and rinse it thoroughly with distilled water;
- 4) Remove excess water using a soft paper towel;
- 5) Place the electrode in a reservoir with potassium chloride (KCl – 3M) solution and let it rest for at least 1 hour;
- 6) After resting, calibrate the instrument.

#### **NOTES:**

- *NEVER keep the electrode in the cleaning solution for more than 30 minutes;*
- *NEVER use abrasives (brushes sandpaper, rough paper) to clean the electrode;*
- *AKSO offers ready-to-use electrode cleaning solutions. Contact us for more information.*

## 8 - MAINTENANCE

---

### DISPOSAL OF BATTERIES AND ELECTRONICS

This product contains batteries and electronic components. Do not dispose of them with other common household waste. Deliver them to the appropriate collection center in accordance with local guidelines.

**Important:** *the correct disposal of electronics and batteries prevents negative effects on the environment and for human health!*

*For more information about disposal services and locations, contact your local municipal office.*

Akso guarantees its instruments against manufacturing defects with the following coverage: 2 years for bench meters, portable and pocket testers and 6 months for electrodes/sensors (unless otherwise indicated).

The warranty period starts from the original date of purchase by the customer. The warranty is only valid if the product is used under normal conditions and in accordance with its operating limits and recommended maintenance procedures.

WARRANTY

2

YEARS


**This instrument comes with a 2-year\* warranty against manufacturing defects and a 6-month\* warranty for the sensor/electrode/probe.**

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