
AKSO[®]

quality you can measure

INSTRUCTION MANUAL

AK150 Max

BENCH PH/EC METER

CONTENTS

1 - SPECIFICATIONS	4
2 - ACCESSORIES	6
3 - OVERVIEW	7
FRONT VIEW	7
VISOR LCD	8
4 - INSTRUCTION OF USE	9
POWER ON/POWER OFF	9
MODE SELECTION AND PARAMETERS (Navigation)	9
MEASUREMENT - pH	10
MEASUREMENT - ORP	10
MEASUREMENT - Conductivity/Resistivity/TDS/Salinity	11
5 - CALIBRATION	13
ADJUSTMENT - pH	13
ADJUSTMENT - ORP	14
ADJUSTMENT - EC	14
6 - ADDITIONAL FUNCTIONS	17
DISPLAY LIGHTING	17
SET CONDUCTIVITY CELL CONSTANT	17
TDS FACTOR ADJUSTMENT	18
FACTORY RESET	18
AUTOMATIC FREEZING OF READING (Auto Hold)	19
7 - SETTINGS	20
SETUP MODE	20
P01 Con - Set conductivity cell constant (K)	20
P02 CAL - Visualization of the specific constant	20
P03 SLP - Set TDS factor	20
P04 buF - Select pH calibration standard	20
P05 SLP - View the slope values from the last calibration	20
P06 AdJ - Set the offset adjustment for temperature measurement	20
P07 Und - Select the temperature measurement unit	20
P08 rES - Setting the numeric display format	21
P09 buF - Select EC adjustment pattern	21
P10 Aut - Enable/disable Automatic reading freezing	21
P11 bL - Enable/Disable backlight	21
P12 Clr - Inactive	21
P13 rSt - Reset	21
8 - MAINTENANCE	22
ELECTRODE CLEANING - pH/ORP	22
DISPOSAL OF BATTERIES AND ELETRONICS	23

1 - SPECIFICATIONS

EC	Measurement range:	0.0 to 19.99 $\mu\text{S/cm}$ 20.0 to 199.9 $\mu\text{S/cm}$ 200 to 1999 $\mu\text{S/cm}$ 2.00 to 19.99 mS/cm
	Resolution:	0.01 $\mu\text{S/cm}$ 0.1 $\mu\text{S/cm}$ 1 $\mu\text{S/cm}$ 0.01 mS/cm
	Accuracy:	$\pm (2\%FS + 1 \text{ digit})$
TDS	Measurement Range:	0.00 to 19.99 ppm 20.0 to 199.9 ppm 200 to 1999 ppm 2.00 to 15.99 ppt
	Resolution:	0.01 ppm 0.1 ppm 1 ppm 0.01 ppt
	Accuracy:	$\pm (2\%FS + 1 \text{ digit})$
Salinity	Measurement Range:	0.00 to 12.00 ppt
	Resolution:	0.01 ppt
	Accuracy:	$\pm (2\%FS + 1 \text{ digit})$
Resistivity	Measurement Range:	0.0 to 100.0 $\text{M}\Omega$
	Resolution:	0.1 $\text{M}\Omega$
	Accuracy:	$\pm (2\%FS + 1 \text{ digit})$
pH	Measurement range:	- 2.00 to 16.00 pH
	Resolution:	0.01 pH
	Accuracy for meter:	$\pm 0.01 \text{ pH}$
ORP	Measurement range:	-1999 to 1999 mV
	Resolution:	1mV
	Accuracy:	$\pm 3 \text{ mV}$

1 - SPECIFICATIONS

Temperature	Measurement range (Meter) :	-5 to 100 °C
	Measurement range (Electrode) :	0 to 80°C
	Measurement range (Cell) :	0 to 80°C
	Resolution:	0.1°C
	Accuracy:	± 0.5°C
Temperature compensation:	Automatic	0 to 80 °C
	Manual:	0 to 80 °C
Temperature compensation (Meter)	Automatic:	-5 to 100°C
	Manual:	-5 to 100°C
Operating temperature of the conductivity cell:		0 to 80°C
Electrode operating temperature:		0 to 80°C
Calibration (Ph):	Calibration (Ph):	4.01, 7.00 and 10.01 pH
	NIST standards:	4 00 6 86 and 9 18 pH
Constant K	Configurable:	0.05 to 0.20
		0.80 to 1.20
		9.80 to 10.20
TDS factor	Configurable:	0.40 to 0.80
Backlight:		On/Off
Reading freeze (HOLD):		Manual
Power supply:		220VAC/12VDC
Operating humidity:		10 to 90 %UR (non-condensing)
Level of protection:		IP65
Dimensions (WxHxD):		216 x 190 x 59 mm
Weight:		950g

Additional functions:

- Automatic calculation of the specific K constant;
- Automatic factory calibration

2 - ACCESSORIES

AK150 Max

Items included with the instrument:

- 1 storage and transport case
- 1 instruction manual
- 1 power supply 220VAC/12VDC – AK1075
- 1 conductivity cell K=1.0 - AK1021-S
- 1 pH Electrode - AK1004-S
- 1 Temperature probe - AK1040
- 1 sample of 1413 $\mu\text{S}/\text{cm}$
- 1 sample of 12.88 mS/cm
- 1 sample of pH 4 solution
- 1 sample of pH 7 solution
- 1 sample of pH 10 solution
- 1 sample of KCl solution


Items sold separately:

- Calibration solution 1413 $\mu\text{S}/\text{cm}$ 250mL – AK4523
- Calibration solution 12.88 mS/cm 250mL – AK4524
- Calibration solution pH 7.00 250mL – AK4507
- Calibration solution pH 4.01 250mL – AK4504
- Calibration solution pH 10.01 250mL – AK4510
- 3M KCl Storage Solution (250 mL bottle) – AK4530
- ORP electrode BNC connector – AK1010

Before use, examine the instrument and accompanying items carefully. If you detect any abnormalities, contact AKSO.

3 - OVERVIEW

FRONT VIEW


AK150 Max:

- 1 – Display
- 2 – On/Off/Set button
- 3 – Mode button
- 4 – Up button
- 5 – Down button
- 6 – CAL/ESC button
- 7 – HOLD/ENTER button

3 - OVERVIEW

DISPLAY LCD

AK150 Max:

- | | |
|--|--|
| 1 - slope display indication | 21 - Manual temperature compensation |
| 2 - pH measurement unit | 22 - Automatic temperature compensation |
| 3 - ORP measurement unit | 23 - Settings menu description |
| 4 - Active reading freeze enabled | 24 - Inactive |
| 5 - Active adjustment/calibration mode | 25 - Inactive |
| 6 - Slope calculation unit | 26 - Salinity measurement mode |
| 7 - Calibrated alkaline adjustment point | 27 - Measurement value |
| 8 - Calibrated neutral adjustment point | 28 - Conductivity measurement stabilization notification |
| 9 - Calibrated acid adjustment point | 29 - Cell constant (K) indication |
| 10 - MΩ measurement unit | 30 - PH/ORP measurement |
| 11 - Parts per million measurement unit | 31 - PH/ORP Measurement stabilization notification |
| 12 - Parts per thousand measurement unit | |
| 13 - uS/cm measurement unit | |
| 14 - mS/cm measurement unit | |
| 15 - Conductivity measurement mode | |
| 16 - Resistivity measurement value | |
| 17 - TDS measurement mode | |
| 18 - Temperature measurement unit in °C | |
| 19 - Temperature measurement unit in °F | |
| 20 - Temperature measurement value | |

4 – INSTRUCTIONS OF USE

POWER ON – POWER OFF

Always connect the power supply to the instrument first. Then, connect the power supply to the outlet.

IMPORTANT NOTICE: The power supply is NOT dual-voltage and operates only at 220 V. Make sure the outlet voltage is correct to avoid damage.

- To turn the instrument on or off, press the button  ;

MODE SELECTION AND PARAMETERS (Navigation)

The instrument features two main measurement groups:

→ Group 1: pH and ORP;

→ Group 2: Salinity, TDS, Resistivity, and Conductivity (EC)

To switch between groups:

Press and hold the  key. The display will alternate between pH/ORP mode and Conductivity mode. Release the key when the desired group is flashing on the display.

To select the parameter (within the group):

With the group already selected (flashing), briefly press the  key (do not hold) to cycle through the options:

- In the pH group: switch between pH and ORP;

- In the Conductivity group: cycle sequentially through Salinity, TDS, Resistivity, and conductivity (EC).

4 – INSTRUCTIONS OF USE

Measurement - pH

- 1) Connect the pH electrode and temperature sensor/probe to the instrument;
- 2) Remove the KCl bottle from the electrode tip, turning it counterclockwise;
- 3) Turn on the instrument by pressing the button  ;
- 4) Rinse the electrode (and temperature probe) with distilled/deionized water and remove excess water using a soft paper towel;
- 5) Immerse the electrode (and temperature probe) in the sample to be analyzed;
- 6) Wait for the reading to stabilize. The display will show 😊;
- 7) After the reading stabilizes, check the pH and temperature values;
- 8) After completing the measurements, rinse the electrode (and temperature probe) with distilled water and store the electrode in the storage bottle with potassium chloride solution (KCl – 3M).

NOTE:

- If a temperature sensor/probe for automatic temperature compensation (ATC) is not connected, the instrument will perform manual temperature compensation (MTC), which can be modified using the button  and  .


Measurement - ORP

- 1) Connect the ORP electrode to the instrument;
- 2) Remove the KCl bottle from the electrode tip;
- 3) Turn on the instrument by pressing the button  ;
- 4) Select the ORP measurement mode by pressing the button  . The display will show the **mV** indicator;
- 5) Rinse the electrode with distilled/deionized water and remove excess water using a soft paper towel;
- 6) Immerse the electrode in the sample to be analyzed;
- 7) Wait for the reading to stabilize. The display will show 😊;
- 8) After the reading stabilizes, check the ORP value;
- 9) After completing the measurements, rinse the electrode with distilled water and store the electrode in the storage bottle with potassium chloride solution (KCl – 3M).

4 – INSTRUCTIONS OF USE

Measurement – Conductivity / Resistivity / TDS / Salinity

- 1) Connect the conductivity cell to the instrument;
- 2) Remove the protection from the cell tip, turning it counterclockwise if it is a cup, or pulling it off if it is a rubber cap.

- 3) Turn on the product with  button;
- 4) Rinse the cell in distilled water and remove excess water using a soft paper towel.
- 5) Immerse the cell in the sample to be analyzed and homogenize it.
- 6) Using the  button, select the desired measurement mode. For each press, the display will show:
 - **EC:** electrical conductivity
 - **RES:** resistivity
 - **TDS:** total dissolved solids
 - **SAL:** salinity

4 – INSTRUCTIONS OF USE

- 7) Wait for the reading to stabilize. The display will show 😊;
- 8) After the reading to stabilize, check the measurement value;
- 9) After completing the measurements, rinse the cell with distilled water and store preferably dry.

NOTES:

- *If a temperature sensor for automatic temperature compensation (ATC) is not connected, the instrument will perform manual temperature compensation (MTC), which can be modified using the button  and  ;*
- *The instrument has the auto-calibration function whenever it is manually restarted after restoring the factory settings. Therefore, calibration with standard solutions should only be performed if a deviation beyond accuracy is detected.*

5 - CALIBRATION

Adjustment - pH

Adjust the instrument's pH measurement upon receiving it and at least once a week. For greater precision, always do this before using the instrument. The instrument can be adjusted to 1, 2 or 3 pH points, respecting the following sequence combinations:

pH7 | pH7 → pH4 | pH7 → pH10 | pH7 → pH4 → pH10

- 1) Connect the pH electrode and temperature sensor/probe to the instrument;
- 2) Remove the KCl bottle from the electrode tip, turning it counterclockwise;
- 3) Turn on the instrument by pressing the button  ;
- 4) Ensure that the measurement selected on the display is pH by viewing the pH unit icon;
- 5) Press the button  to access the pH measurement adjustment mode;
- 6) The display will show an indicator  at the top of the screen, signaling that the meter is in calibration mode;
- 7) Rinse the electrode (and temperature probe) in distilled/deionized water and remove excess water using a soft paper towel;
- 8) Using a separate bottle, take a sufficient portion of the pH solutions, enough to cover the electrode tip. **NEVER** insert the electrode into the original vials of 250mL.
- 9) Immerse the electrode (and temperature probe) in the pH buffer solution, agitating it gently to homogenize it;
- 10) Wait for the reading to stabilize. The notification  will appear on the screen;
- 11) After the reading stabilizes, press the button  to save the adjustment;
 - To calibrate another pH point, repeat steps 7 to 10 using the corresponding buffer solution;
 - To exit adjustment mode, press the button . The instrument will return to measurement mode.
- 12) After completing the adjustment, the instrument will display the pH points that have been adjusted next to the pH value;
- 13) After finishing using, rinse the electrode (and temperature probe) with distilled/deionized water and store the electrode in the storage bottle with potassium chloride solution (KCl – 3M).

5 - CALIBRATION

NOTE:

- Always begin the adjustment at the pH 7.00 point and, after completing the procedure, discard the used portions of the solutions;
- When saving the adjustment at point pH 4.01 or pH 10.01, the display will show the electrode slope percentage. The indicator **Slope** will appear;
- The slope indicates the electrode's performance and should be between 80% and 115%. If the slope is below 80% or above 115% clean the electrode, refer to **7 - MAINTENANCE>ELECTRODE CLEANING – pH/ORP**.
- If the message **Err** appears on the display when trying to save an adjustment, check the pH solutions and the electrode and repeat the procedure;
- When adjusting the 3 pH points, the instrument will automatically return to measurement mode.

Adjustment - ORP

The ORP electrode does not require adjustment. However, you can check the electrode's performance and the accuracy of the readings using an ORP standard solution.

Adjustment - EC

The calibration with the standard solution should be performed only with the solution closest to the operational measurement range in which conductivity is being measured.

The instrument has two adjustment buffers that can be selected through the **P09 bUF** menu: Std (Standard) and Cst (Custom).

Std (Standard) – Fixed set points at values: 84 $\mu\text{S}/\text{cm}$, 1413 $\mu\text{S}/\text{cm}$, 5 mS/cm , and 12.88 mS/cm .

CSt (Custom) – Customizable set points using the arrows after entering calibration mode.

To perform the adjustment in **Std** mode, follow the instructions below:

- 1) Connect the conductivity cell to the instrument.
- 2) Remove the protective bottle/cap from the cell tip.
- 3) Turn on the instrument by pressing the button .
- 4) Ensure that the measurement selected on the display is electrical conductivity, verifying the **EC** icon on the display (if the instrument is set to another parameter, entering adjustment mode will automatically display the conductivity measurement value).
- 5) Press the  button to access the adjustment mode.

5 - CALIBRATION

- 6) The indication **Cal** will appear at the top of the screen, signaling that the meter is in calibration mode.
- 7) Rinse the conductivity cell with distilled water and remove any excess water using a soft paper towel.
- 8) Pour into a separate container a portion of the solution to be used (84 $\mu\text{S}/\text{cm}$, 1413 $\mu\text{S}/\text{cm}$, 5 mS/cm , or 12.88 mS/cm). Make sure there is enough volume to cover the cell tip. **NEVER** place the electrode into the original solution bottles.
- 9) Immerse the conductivity cell in the standard solution, gently stirring it to homogenize.
- 10) Wait for the reading to stabilize. Then, using the  button, confirm the adjustment. The instrument will automatically identify the appropriate solution and perform the calibration.

NOTE:

- If the value is outside the expected adjustment range for the nearest point, the display will show the message "Err", indicating that there is a problem (contaminated buffer, defective conductivity cell, or some error in the procedure).
- *The adjustment of the conductivity parameter affects the other calculations performed by the meter.*

5 - CALIBRATION

To make adjustments in **Cst** mode, follow the guidelines below:

- 1) In measurement mode, press and hold the  button to enter the setup mode. The display will show the indication P01;
- 2) Using the buttons  and , navigate to the parameter **P09 bUF** and briefly press the  button to access it; Navigate to parameter **P09 bUF**. It is set to Std by default; change it to CSt.
- 3) Turn off the instrument by pressing the  button
- 4) Connect the conductivity cell to the instrument and next remove the protective cap from the tip of the cell
- 5) Turn on the instrument by pressing the  button;
- 6) Ensure that the measurement selected on the display is electrical conductivity by checking the EC icon on the display (if the instrument is set to another parameter, the conductivity measurement value will automatically be displayed when entering the adjustment mode);
- 7) Press the  button to enter the adjustment mode;
- 8) The indicator  will appear at the top of the screen signaling that the meter is in calibration mode. Note that the displayed value refers to the conductivity measurement without any adjustments;
- 9) Rinse the conductivity cell with distilled water and remove the excess water using a soft paper towel;
- 10) Prepare a portion of the solution to be used in a separate container. Ensure there is enough volume to cover the tip of the cell. NEVER insert the electrode into the original 250mL containers.
- 11) Immerse the conductivity cell in the standard solution, gently stirring to homogenize it;
- 12) Wait for the reading to stabilize. Then, using the buttons  and , adjust the measured value to the actual value of your solution. If necessary, press and hold the buttons to correct the value more quickly.
- 13) Once the reading has stabilized, press the  button to save the adjustment. The instrument will automatically return to measurement mode with the corrected value.

6 - ADDITIONAL FUNCTIONS


DISPLAY LIGHTING

To activate/deactivate the display lighting:

- 1) In measurement mode, press and hold the  button to enter the setup mode.
The display will show the indication P01;
- 2) Using the buttons  and  navigate to the parameter **P11 bL** and briefly press the  button to access it;
- 3) Using the buttons  and  select the option **On** to activate or **OFF** to deactivate. Press the  button to confirm.
- 4) Press the  button to return to measurement mode.

SET CONDUCTIVITY CELL CONSTANT

The exact value of the cell constant is usually displayed on the cable, as shown in the figure:

- 1) In measurement mode, press and hold the  button to enter the setup mode.
The display will show the indication P01;
- 2) Using the buttons  and  navigate to the parameter **P01 Con** and briefly press the  button to access it;
- 3) Using the buttons  and  , increment or decrement the value according to your conductivity cell. Then, press the  button to confirm.
- 4) The product will return to configuration mode. Check if the option saved in memory is desired one.

6 - ADDITIONAL FUNCTIONS

TDS FACTOR ADJUSTMENT

The AK150 has the function of estimating the total dissolved solids in the solution by converting conductivity using a known constant. The TDS factor depends on the sample's reference salt and can be adjusted as needed by following the instructions below:

- 1) In measurement mode, press and hold the  button to enter the setup mode.
The display will show the indication P01;
- 2) Using the buttons  and , navigate to the parameter **P03 tdS** and briefly press the  button to access it;
- 3) Using the buttons  and , increment or decrement the value according to your reference salt. Then, press the  button to confirm.
- 4) The product will return to configuration mode. Check if the option saved in memory is the desired one.
- 5) Press the  button to return to measurement mode.

FACTORY RESET

Whenever an adjustment is made incorrectly or a complete restoration of the equipment's firmware is required, follow the instructions below:

- 1) In measurement mode, press and hold the  button to access the configuration mode. The display will show the indication P01 ;
- 2) Using the buttons  and , navigate to the parameter **P13 rSt** and briefly press the button  to access it;
- 3) "Yes" will be displayed on the screen. To confirm the reset, press the  button.
- 4) After the product automatically returns to the measurement screen, turn the instrument off and then on again.

6 - ADDITIONAL FUNCTIONS

AUTOMATIC FREEZING OF READING (Auto-Hold)

When enabled, the meter will automatically set the reading value on the device's display after the measurement has stabilized. To enable or disable Auto-Hold:


- 1) While in measurement mode, press and hold the button  to enter configuration mode. The display will show P01;
- 2) Using the buttons  and  navigate to the parameter **P10 Aut** and briefly press the button  to access it;
- 3) Using the buttons  and  select the option **On** to activate or **OFF** to deactivate. Press the button  to confirm.
- 4) The product will return to the setup mode. Verify that the option saved in memory is the desired one.
- 5) Press the button  to return to measurement mode.

7 - SETTINGS

SETUP MODE

- To access the instrument's configurable parameter menu, in measurement mode, press and hold the button . The display will show the indication P01;
- To navigate between parameters, use the buttons  e .
- To access a parameter, press the button .
- To toggle between the parameter configuration options use the buttons  and .
- To confirm a setting, press the button . The instrument will return to the parameter selection;
- To go back to a previous level, press the button .

The figure below shows how the parameters are displayed:

Available parameters for configuration:

P01 Con – Set conductivity cell constant (K)

→ **0.10** (0.05 to 0.20)

→ **1.00** (0.80 to 1.20)

→ **10.00** (9.80 to 10.20)

P02 CAL – Visualization of the specific constant

This menu is dedicated to show the constant that the instrument is using in its measurements after user calibration.

P03 tds – Set TDS factor

→ **0.30 to 1.00**

P04 buF – Select pH calibration standard

→ **USA** (4.01, 7.00 and 10.01 pH)

→ **NST** (4.01, 6.86 and 9.18 pH)

P05 SLP – View the slope values from the last calibration

P06 AdJ – Set the offset adjustment for temperature measurement

→ **± 5.0°C** (from the measured value)

P07 Und – Select the temperature measurement unit

→ **°C** (Celsius)

→ **°F** (Fahrenheit)

7 - SETTINGS

P08 rES – Setting the numerical display format

P09 buF – Select EC adjustment pattern

→ **Std** (standart) 84 $\mu\text{S/cm}$, 1413 $\mu\text{S/cm}$, 5 mS/cm , 12.88 mS/cm

→ **CSt** (custom) customize

P10 Aut – Enable/disable Automatic reading freeze

→ **On** (enable)

→ **OFF** (disable)

P11 bL – Enable/disable backlight

→ **On** (enable)

→ **OFF** (disable)

P12 CLr – Inactive

→

P13 rSt – Restore factory settings

→ **YES** (perform the reset)

Factory Settings of the Instrument

P01 Con	- Set Conductivity Cell Constant	1.00
P02 CAL	- Visualization of the Specific Constant	***
P03 tdS	- Set TDS factor	0.80
P04 buF	- pH Calibration Standard	USA
P05 SLP	- View Slope Values	***
P06 AdJ	- Set Temperature Offset Adjustment	0.0
P07 Und	- Select the temperature measurement unit	$^{\circ}\text{C}$
P08 rES	- Setting the numerical display format	0.01
P09 buF	- Select EC adjustment pattern	Std
P10 Aut	- Enable/disable automatic reading freeze	***
P11 bL	- Enable/Disable backlight	ON
P12 CLr	- Inactive	***
P13 rSt	- Reset	***

Ensure that the temperature sensor/probe is connected to perform the temperature measurement offset adjustment.

8 - MAINTENANCE

ELECTRODE CLEANING - pH / ORP

To ensure measurement accuracy and prolong the life of the pH/ORP electrode, perform regular cleaning (at least every two weeks) or when the electrode slope falls below 90%:

- 1) Prepare a container with an electrode cleaning solution (5% pepsin in HCl 0.1M solution);
- 2) Immerse the electrode tip in the solution for 20 to 25 minutes;
- 3) After this period, remove the electrode from the solution and rinse it thoroughly with distilled water;
- 4) Remove excess water using a soft paper towel;
- 5) Place the electrode in a reservoir with potassium chloride (KCl – 3M) solution and let it rest for at least 1 hour;
- 6) After resting, calibrate the instrument.

NOTES:

- *NEVER keep the electrode in the cleaning solution for more than 30 minutes;*
- *NEVER use abrasives (brushes sandpaper, rough paper) to clean the electrode;*
- *AKSO offers ready-to-use electrode cleaning solutions. Contact us for more information.*

8 - MAINTENANCE

DISPOSAL OF BATTERIES AND ELECTRONICS

This product contains batteries and electronic components. Do not dispose of them with other common household waste. Deliver them to the appropriate collection center in accordance with local guidelines.

Important: *the correct disposal of electronics and batteries prevents negative effects on the environment and for human health!*

For more information about disposal services and locations, contact your local municipal office.

Akso guarantees its instruments against manufacturing defects with the following coverage: 2 years for bench meters, portable and pocket testers and 6 months for electrodes/sensors (unless otherwise indicated).

The warranty period starts from the original date of purchase by the customer. The warranty is only valid if the product is used under normal conditions and in accordance with its operating limits and recommended maintenance procedures.

WARRANTY

2

YEARS

This instrument comes with a 2-year* warranty against manufacturing defects and a 6-month* warranty for the sensor/electrode/probe.

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