

## PH PROBE

# AK1001 PH ELECTRODE FOR LIQUIDS

The AK1001 pH electrode is a combined electrode, ideal for a wide range of general applications in aqueous solutions. It has an integrated temperature sensor, which allows automatic temperature compensation, facilitating handling and providing faster and more accurate measurements.



## SPECIFICATIONS

Measuring range:	0 to 14 pH
Operating temperature:	0 to 80 °C
Junction:	Fiber
Internal electrolyte:	Ag / AgCl
Connection:	BNC
Material:	PC (polycarbonate)
Dimensions:	Electrode: 160 x 13 mm (Ø) Cable: 1m
Weight	60g

**6** months warranty  
(PROBE)