

## PH PROBE

# AK1004 PH ELECTRODE FOR LIQUIDS

The AK1004 pH electrode is a combined electrode, ideal for a wide range of general applications in aqueous solutions. It has an integrated temperature sensor, which allows automatic temperature compensation, facilitating handling and providing faster and more accurate measurements.



## SPECIFICATIONS

|                        |   |
|------------------------|---|
| Measuring range:       | 0 to 14 pH                              |
| Operating temperature: | 0 to 90 °C                              |
| Junction:              | Ceramic                                 |
| Internal electrolyte:  | Ag / AgCl                               |
| Connection:            | BNC                                     |
| Material:              | Glass                                   |
| Dimensions:            | Electrode: 155 x 12 mm (Ø)<br>Cable: 1m |
| Weight                 | 65g                                     |

**6** months warranty  
(PROBE)